

Managing skin tears in practice

It is estimated that prevalence of skin tears may be underreported and in fact be greater than pressure ulcers¹ – to date, no prevalence data is available for the UK, so the cost to patients and the NHS is not fully known²

What are skin tears?

- Skin tears are acute wounds caused by shear, friction or trauma, resulting in separation of the skin layers³
- Skin tears can be full or partial thickness and can occur anywhere on the body – most commonly seen on the hands, arms and lower legs
- 70–80% of skin tears occur on hands or arms⁴
- It is estimated that prevalence of skin tears may be underreported and in fact be greater than pressure ulcers¹ – to date, no prevalence data is available for the UK, so the cost to patients and the NHS is not fully known²
- A US study reported 1.5 million skin tears affect in-patients every year⁴
- The ageing population means that incidence of skin tears is increasing (elderly patients have fragile skin and are at increased risk)⁴
- Skin must be protected in at-risk patients and skin tears managed to avoid further damage and complication⁴
- Skin tears can be painful and distressing for the patient⁴

Skin tear risk assessment (patient, wound, environment)⁵

Risk categories

- **Skin:** extremes of age, dry/fragile skin, previous skin tear
 - **Mobility:** history of fall, impaired mobility, dependent activities of daily living, mechanical trauma
 - **General health:** comorbidities, polypharmacy, impaired cognition (sensory, visual, auditory), malnutrition

At Risk

If patient has **any** identified risk factors


NO

Reassess when patient's condition changes

YES

Implement risk reduction programme checklist and skincare regimen

Managing skin tears step by step (adapted from Wounds UK⁵)



Skin tear categories

Type 1
No skin loss



Type 2
Partial flap loss



Type 3
Total flap loss






Dressing selection for managing skin tears

Dressing selection is a key element of managing skin tears and it is important to select the appropriate dressing with treatment goals in mind. As such, the ideal dressing for managing skin tears should:¹

- Control bleeding
- Be easy to apply
- Provide a protective anti-shear barrier
- Optimise the physiological healing environment (e.g. moisture, bacterial balance, temperature, pH)
- Be flexible and mould to contours
- Provide secure, but not aggressive, retention
- Afford extended wear time
- Not cause trauma on removal
- Optimise quality of life and cosmesis
- Be cost-effective

Dressing selection for managing skin tears

STAR skin tear classification system	ISTAP skin tear classification system	Skin tear treatment options as recommended by ISTAP	Acelyty™ dressing options
<p>Category 1A and 1B</p> <p>1A: A skin tear where the edges can be realigned to the normal anatomical position (without undue stretching) and the skin or flap colour is not pale or darkened</p> <p>1B: A skin tear where the edges can be realigned to the normal anatomical position (without undue stretching) and the skin or flap colour is pale or darkened. (Image represents Star 1B)</p>	<p>Type I: Skin tear without tissue loss No skin loss; linear or flap tear, which can be repositioned to cover the wound bed</p> 	<p>Based on assessment Control bleeding; approximate edges. Cover wound with a silicone contact layer. Apply appropriate secondary dressing when required, such as a non-adhesive or silicone foam, depending on wound exudate and location.</p>	<p>ADAPTIC TOUCH™ Non-Adhering Silicone Dressing KERRAFOAM™ Gentle Border with EXUSAFE™ Technology KERRALITE COOL™ Dressing or KERRALITE COOL™ Border Dressing for painful nil to low exuding wounds.</p>
<p>Category 2A and 2B</p> <p>2A: A skin tear where the edges cannot be realigned to the normal anatomical position and the skin or flap is not pale or darkened</p> <p>2B: A skin tear where the edges cannot be realigned to the normal anatomical position and the skin or flap colour is pale or darkened. (Image represents Star 2B)</p>	<p>Type II: Partial flap loss Flap cannot be repositioned to cover the wound</p> 	<p>Control bleeding; approximate edges. Cover wound with a silicone contact layer. Apply appropriate secondary dressing when required, such as a non-adhesive or silicone foam, depending on wound exudate and location.</p>	<p>ADAPTIC TOUCH™ Non-Adhering Silicone Dressing KERRAFOAM™ Gentle Border with EXUSAFE™ Technology KERRALITE COOL™ Dressing or KERRALITE COOL™ Border Dressing for painful nil to low exuding wounds.</p>
<p>Category 3</p> <p>A skin tear where the skin flap is completely absent</p>	<p>Type III: Total flap loss Entire wound bed is exposed</p> 	<p>Control bleeding; cover wound with a non-adhering silicone contact layer. Apply appropriate secondary dressing when required, such as a non-adhesive or silicone foam, depending on wound exudate and location.</p>	<p>ADAPTIC TOUCH™ Non-Adhering Silicone Dressing KERRAFOAM™ Gentle Border with EXUSAFE™ Technology KERRALITE COOL™ Dressing or KERRALITE COOL™ Border Dressing for painful nil to low exuding wounds.</p> <p>For bleeding, partial or total flap loss: PROMOGRAN™ Protease Modulating Matrix or PROMOGRAN PRISMA™ Wound Balancing Matrix when at risk of infection*</p>

* Apply as a primary wound contact layer, then cover with an appropriate secondary dressing. For low or no exudate, use saline to moisten the matrix and initiate transformation into gel. Note: PROMOGRAN™ Matrix and PROMOGRAN PRISMA™ Matrix: If gel has not biodegraded, it is not necessary to remove.

References

1. Stephen Haynes J, Carville K (2011) Skin Tears Made Easy. *Wounds International* 2(4): 1-6.
2. Bianchi J (2012) Preventing, assessing and managing skin tears. *Nursing Times* 108: 13, 12-6.
3. LeBlanc K et al (2013) International Skin Tear Advisory Panel: a tool kit to aid in the prevention, assessment, and treatment of skin tears using a Simplified Classification System ©. *Adv Skin & Wound Care* 26(10): 451.
4. Hebert G (2016) No more skin tears. *Wound Care Advisor* 5:2.
5. Wounds UK (2015) All Wales guidance for the prevention and management of skin tears. Available at: www.welshwoundnetwork.org/files/8314/4403/4358/content_11623.pdf.

Note: Specific indications, contraindications, warnings, precautions and safety information may exist for Systagenix and KCI (Acelyty companies) products. Please consult a healthcare provider and product instructions for use prior to application.
Copyright 2019 KCI Licensing, Inc. All rights reserved. Unless otherwise noted, all trademarks designated herein are proprietary to KCI Licensing, Systagenix Wound Management, Ltd., or Crawford Healthcare, Ltd. PRA-PM-EU-00178 (05/19).